



PRODUCT DATA SHEET

Product: Fas Ligand, soluble Enhancer

Cat. No.: FL-104 (50 µg)

Specificity:

Enhances biological activity of TRAIL, soluble (human) (recombinant) and FasL, soluble (human) (recombinant).

Format:

Liquid. In PBS, pH 7.4, containing 50% glycerol. Contains no preservatives.

Storage:

Store at -20°C long term. After opening, prepare appropriate aliquots and store them at -20°C. Avoid repeated freeze/thaw cycles.

Applications and Suggested Dilutions:

- Negative control in experiments with recombinant proteins fused to human IgG1 Fc.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.