


Mouse Isotyping Kit

Catalog Number: KT-789

- For rapid identification of mouse immunoglobulin isotypes, subtypes and light chains.
- Each kit is suitable for running 20 unique samples.
- Conveniently store this product at room temperature.
- Kit includes:

Components
20 red strips to identify IgG1, IgG2a, IgG2b and IgG3 subtypes
20 blue strips to identify kappa light chain, lambda light chain, IgA and IgM
Sample Buffer



For *In Vitro* Research Use Only. Not for use in diagnostic procedures in the U.S.
For additional information please visit K-ASSAY.com

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168
Tel: (206) 575-8068 Fax: (206) 575-8094
Email: LifeScience@k-assay.com
www.k-assay.com