
LPO-CC (Lipid Peroxides) Kit

Catalog Number: CC-004



- No sample extraction required for serum or plasma.
- Specifically quantitates lipid peroxides. (Does not measure metabolites of lipid peroxides.)
- Conveniently store this product at 4°C.
- Manual or automated procedures available.

Kit includes:

Component	Size
Reagent 1: Enzyme Reagent	55 mL x 1
Reagent 2A: Lyophilized Chromogen Reagent	Lyoph. x 1
Reagent 2B: Buffer Solution (liquid)	110 mL x 1
Calibrator (liquid)	10 mL x 1



For *In Vitro* Research Use Only. Not for use in diagnostic procedures in the U.S. For additional information please visit K-ASSAY.com

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168 Tel: (206) 575-8068 Fax: (206) 575-8094 Email: LifeScience@k-assay.com www.k-assay.com