


Mouse Embryo Cryopreservation Kit

Catalog Number: KT-102

- For the cryopreservation and thawing of mouse embryos.
- Kit is suitable for cryopreserving about 8,000 mouse embryos.
- Conveniently store this product at 4°C.
- Kit includes:

Components for Cryopreservation

Component	Size
TYH Medium	5 mL x 1
mW Medium	5 mL x 1
PB1 Solution	5 mL x 1
1M DMSO Solution	2 mL x 1
DAP213 Preservation Solution	2 mL x 1



Components for Thawing

Component	Size
0.25M Sucrose Solution	2 mL x 5
mW Medium	5 mL x 5

For *In Vitro* Research Use Only. Not for use in diagnostic procedures in the U.S.

For additional information please visit K-ASSAY.com

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168
Tel: (206) 575-8068 Fax: (206) 575-8094
Email: LifeScience@k-assay.com
www.k-assay.com