

PRODUCT DATA SHEET

Product: Rabbit Reference Serum

Cat. No.: KT-098 (1.0 mL)

Description:

Rabbit Reference Serum is an assayed control material for use with **K-ASSAY** Rabbit Haptoglobin (KT-352), CRP (KT-097) and Transferrin (KT-351) ELISA kits.

Format:

1 mL liquid rabbit reference serum with 0.1% sodium azide. Assayed values for rabbit CRP, rabbit Haptoglobin, and rabbit Transferrin are available (assayed values were obtained using the corresponding **K-ASSAY** ELISAs). See Control Certificate for the lot specific concentrations.

Procedure:

Serum was collected from rabbits.

Storage:

For storage longer than 7 days, keep frozen. Storage less than 7 days can be at 4°C. Avoid multiple freeze/thaw cycles.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.