



KAMIYA BIOMEDICAL COMPANY

Monkey Cystatin C ELISA

For the quantitative determination of Cystatin C (CYS) in monkey biological fluid

Cat. No. KT-666

For research use only.



PRODUCT INFORMATION

Monkey Cystatin C ELISA Cat. No. KT-666

INTENDED USE

The Monkey CYS ELISA is a highly sensitive two-site enzyme-linked immunoassay (ELISA) for the quantitative determination of CYS in monkey biological fluid. For research use only.

INTRODUCTION

Cystatin C is a small cysteine proteinase inhibitor present in body fluids. Studies have shown Cystatin C levels to be directly correlated with the glomular filtration rate.

PRINCIPLE

The principle of the double antibody sandwich ELISA is represented in Figure 1. In this assay the CYS present in the sample reacts with the anti-CYS antibodies which have been adsorbed to the surface of polystyrene microtiter wells. After the removal of unbound proteins by washing, anti-CYS antibodies conjugated with horseradish peroxidase (HRP) are added. These enzyme-labeled antibodies form complexes with the previously bound CYS. Following another washing step, the complexes are assayed by the addition of a chromogenic substrate, 3,3',5,5'-tetramethylbenzidine (TMB). The quantity of bound enzyme is proportional to the concentration of CYS in the sample tested; thus, the absorbance, at 450 nm, is a measure of the concentration of CYS in the test sample. The quantity of CYS in the test sample can be interpolated from the calibration curve constructed from the calibrators, and corrected for sample dilution.

Figure 1.	Anti-CYS Antibodies Bound To Solid Phase	
	Calibrators and Samples Added	
	CYS * Anti-CYS Complexes Formed	
	Unbound Sample Proteins Removed	
	Detection Antibody Added	
	Complexes Formed	
	Unbound Protein Removed	
	HRP Streptavidin Added	
	Complexes Formed	
	Unbound HRP Removed	
	Chromogenic Substrate Added	
	Determine Bound Enzyme Activity	

COMPONENTS

- 1. Diluent One bottle containing 60 mL of a 1X diluent running buffer.
- 2. Wash Solution Concentrate One bottle containing 50 mL of a 20X concentrated wash solution.

3. Detection Antibody 100X

One vial containing 150 μ L of affinity purified anti-Monkey CYS antibody conjugated with horseradish peroxidase in a stabilizing buffer.

- HRP-Streptavidin 100X
 One vial containing 150 μL of Horseradish Peroxidase conjugated streptavidin in a stabilizing buffer.
- TMB Substrate Solution One vial containing 12 mL of TMB and hydrogen peroxide in citric acid buffer at pH 3.3.
- 6. Stop Solution

One vial containing 12 mL of 0.3 M sulfuric acid. WARNING: Avoid contact with skin.

7. Microtiter Plate

Twelve removable eight-well micro strips in well holder frame. Wells are coated with affinity-purified anti-monkey CYS.

8. Monkey CYS Calibrator

One vial containing a lyophilized Monkey CYS Calibrator. The calibrator used in this kit is of human origin. The use of a human calibrator allows export of this kit.

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes (2 µL to 200 µL) for making and dispensing dilutions
- Test tubes
- Microplate washer/aspirator
- Distilled or de-ionized H₂O
- Microplate reader
- Assorted glassware for the preparation of reagents and buffer solutions
- Timer
- Orbital Shaker

PRECAUTIONS

- 1. Read the instructions carefully before beginning the assay.
- 2. This kit is for research use only.
- 3. Great care has been taken to ensure the quality and reliability of this product. However, it is possible that in certain cases, unusual results may be obtained due to high levels of interfering factors.
- 4. No additives or preservatives are necessary to maintain the integrity of the specimen. Avoid azide contamination.
- 5. Azide and thimerosal at concentrations higher than 0.1% inhibit the enzyme reaction.
- 6. Other precautions:
 - > Do not interchange kit components from different lots.
 - > Do not use kit components beyond the expiration date.
 - Protect reagents from direct sunlight.
 - Do not pipette by mouth.
 - > Do not eat, drink, smoke or apply cosmetics where reagents are used.
 - > Avoid all contact with the reagents by using gloves.
 - Stop solution contains diluted sulfuric acid. Irritation to eyes and skin is possible. Flush with water after contact.

REAGENT PREPARATION

- 1. Diluent Concentrate Ready to use as supplied.
- 2. Wash Solution Concentrate

The Wash Solution supplied is a 20X concentrate and must be diluted 1:20 with distilled or de-ionized water. Crystal formation in the concentrate is not uncommon when storage temperatures are low. Warming of the concentrate to 30-35°C before dilution can dissolve crystals.

3. Detection Antibody 100X

Calculate the required amount of working conjugate solution for each microtiter plate test strip by adding 10 μ L detection antibody to 990 μ L of 1X Diluent for each test strip to be used for testing. Mix uniformly, but gently. Avoid foaming.

4. HRP-Streptavidin 100X

Calculate the required amount of working conjugate solution for each microtiter plate test strip by adding 10 µL HRP-Streptavidin to 990 µL of 1X Diluent for each test strip to be used for testing. Mix uniformly, but gently. Avoid foaming.

- 5. TMB Substrate Solution Ready to use as supplied.
- 6. Stop Solution Ready to use as supplied.
- 7. Microtiter Plate

Ready to use as supplied. Unseal Microtiter Pouch and remove plate from pouch. Remove all strips and wells that will not be used in the assay and place back in pouch and re-seal along with desiccant.

8. Monkey CYS Calibrator

Add 1.0 mL of distilled or de-ionized water to the lyophilized Monkey CYS Calibrator and mix gently until dissolved. It is at a concentration of 1.236 µg/mL (the reconstituted calibrator should be aliquoted and frozen if future use is intended). CYS Calibrators need to be prepared immediately before use (see chart below). Mix well between each step. Avoid foaming.

Calibrator	Concentration (ng/mL)	Calibrator Volume added to 1X Diluent	Volume of 1X Diluent
6	4	2 μL Monkey CYS Calibrator	616 μL
5	2	0.3 mL Calibrator 6	0.3 mL
4	1	0.3 mL Calibrator 5	0.3 mL
3	0.5	0.3 mL Calibrator 4	0.3 mL
2	0.25	0.3 mL Calibrator 3	0.3 mL
1	0.125	0.3 mL Calibrator 2	0.3 mL
0	0		0.6 mL

STORAGE AND STABILITY

1. Complete Kit

The expiration date for the kit is stated on the outer label. The recommended storage temperature is 4°C. Note: See long term storage recommendations below for the Monkey CYS Calibrator.

2. Diluent

The 1X Diluent Concentrate is stable until the expiration date and should be stored at 4°C.

3. Wash Solution

The 20X Wash Solution Concentrate is stable until the expiration date. The 1X working solution is stable for at least one week from the date of preparation. Both solutions can be stored at room temperature (RT, 16-25°C) or at 4°C.

4. Detection Antibody 100X

Undiluted Biotin conjugated anti-CYS should be stored at 4°C and diluted immediately prior to use. The working conjugate solution is stable for up to 1 hour when stored in the dark.

5. HRP-Streptavidin 100X

Undiluted horseradish peroxidase conjugated streptavidin should be stored at 4°C and diluted immediately prior to use. The working conjugate solution is stable for up to 1 hour when stored in the dark.

TMB Substrate Solution The TMB Substrate Solution should be stored at 4°C and is stable until the expiration date.

7. Stop Solution

The Stop Solution should be stored at 4°C and is stable until the expiration date.

8. Microtiter Plate

Anti-monkey CYS coated wells are stable until the expiration date and should be stored at 4°C in the sealed foil pouch with a desiccant pack.

9. Monkey CYS Calibrator

The lyophilized Monkey CYS Calibrator should be stored at 4°C or frozen until reconstituted. The reconstituted calibrator should be aliquoted and stored frozen. Avoid multiple freeze/thaw cycles. The working calibrator solutions should be prepared immediately prior to use.

INDICATIONS OF INSTABILITY

If the test is performing correctly, the results observed with the calibrator solutions should be within 20% of the expected values.

SPECIMEN COLLECTION AND HANDLING

Blood should be collected by venipuncture and the serum separated from the cells, after clot formation, by centrifugation. For plasma samples, blood should be collected into a container with an anticoagulant and then centrifuged. Care should be taken to minimize hemolysis, excessive hemolysis can impact your results. Assay immediately or aliquot and store samples at -20°C. Avoid repeated freezing/thawing.

For any sample that might contain pathogens, care must be taken to prevent contact with open wounds. No additives or preservatives are necessary to maintain the integrity of the specimen. Avoid azide contamination.

ASSAY PROTOCOL

Dilution of Samples

The assay for quantification of CYS in samples requires that each test sample be diluted before use. For a single step determination a dilution of 1:1,000 is appropriate for most serum/plasma samples. For absolute quantification of samples that yield results outside the range of the calibration curve, a lesser or greater dilution might be required. If unsure of sample level, a serial dilution with one or two representative samples before running the entire plate is highly recommended.

To prepare a 1:1,000 dilution of sample, transfer 2 μ L of sample to 198 μ L of 1X Diluent. This gives you a 1:100 dilution. Next, dilute the 1:100 by transferring 30 μ L of sample to 270 μ L of 1X diluent. This gives you a 1:1,000 dilution. Mix thoroughly at each stage. Dilute samples immediately prior to use.

Procedure

- 1. Bring all reagents to RT before use.
- 2. Pipette 100 µL of
 - Calibrator 0 (0.0 ng/mL) in duplicate Calibrator 1 (0.125 ng/mL) in duplicate Calibrator 2 (0.25 ng/mL) in duplicate Calibrator 3 (0.5 ng/mL) in duplicate Calibrator 4 (1 ng/mL) in duplicate Calibrator 5 (2 ng/mL) in duplicate Calibrator 6 (4 ng/mL) in duplicate
- 3. Pipette 100 µL of diluted sample (in duplicate) into pre-designated wells.
- 4. Incubate the Microtiter Plate at 22°C (RT) for sixty (60 ± 2) minutes while shaking on an orbital shaker. Keep plate covered and level during incubation.
- 5. Following incubation, aspirate the contents of the wells.
- 6. Completely fill each well with appropriately diluted Wash Solution and aspirate. Repeat three times, for a total of four washes. If washing manually; completely fill wells with wash buffer, invert the plate and pour/shake out the contents in a waste container. Follow this by sharply striking the wells on absorbent paper to remove residual buffer. Repeat 3 times for a total of 4 washes.

- Pipette 100 μL of appropriately diluted detection antibody to each well. Incubate at 22°C (RT) for twenty (20 ± 2) minutes while shaking on an orbital shaker. Keep plate covered in the dark and level during incubation.
- 8. Wash and blot the wells as described in Steps 5 and 6.
- 9. Pipette 100 µL of appropriately diluted HRP-streptavidin to each well. Incubate at 22°C (RT) for twenty (20 ± 2) minutes while shaking on an orbital shaker. Keep plate covered in the dark and level during incubation.
- 10. Wash and blot the wells as described in Steps 5 and 6.
- 11. Pipette 100 µL of TMB Substrate Solution into each well.
- 12. Incubate in the dark at RT for precisely ten (10) minutes while shaking on an orbital shaker.
- 13. After ten (10) minutes, add 100 µL of Stop Solution to each well.
- 14. Determine the absorbance at 450 nm of the contents of each well. Zero the plate reader to air.

The absorbance of the final reaction mixture can be measured up to 2 hours after the addition of the Stop Solution. However, good laboratory practice dictates that the measurement be made as soon as possible.

RESULTS

- 1. Subtract the average background value from the test values for each sample.
- 2. Using the results observed for the calibrators construct a calibration curve. The appropriate curve fit is that of a fourparameter logistics curve. A second order polynomial (quadratic) or other curve fits may also be used.
- 3. Interpolate test sample values from calibration curve. Correct for sample dilution factor to arrive at CYS concentration in original sample.

QUALITY CONTROL

In accord with good laboratory practice, the assays for specific CYS require meticulous quality control. Each laboratory should use routine quality control procedures to establish inter- and intra-assay precision and performance characteristics.

LIMITATION OF THE PROCEDURE

- Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the information contained in the package insert instructions and with adherence to good laboratory practice.
- 2. Factors that might affect the performance of the assay include proper instrument function, cleanliness of glassware, quality of distilled or de-ionized water, and accuracy of reagent and sample pipettings, washing technique, incubation time or temperature.

FOR RESEARCH USE ONLY

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168 Tel: (206) 575-8068 Fax: (206) 575-8094 Email: LifeScience@k-assay.com www.k-assay.com