



PRODUCT DATA SHEET

Product: Anti-Aggrecan mAb, clone 8E10

Cat. No.: MC-1000 (100 µL)

Specificity:

Reacts with the GAG linkage region of aggrecan.

Species Reactivity:

Human, bovine, and shark. Other species not tested.

Ig Isotype:

Mouse IgG2a

Immunogen:

Human articular cartilage aggrecan

Format:

100 µL of ascites.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- ELISA
- Immunohistochemistry
- Immunoprecipitation
- Western blot

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.