

PRODUCT DATA SHEET

Product: Anti-Aggrecan mAb, clone 7E1

Cat. No.: MC-1001 (100 µL)

Specificity: Reacts with the G3 domain of aggrecan.

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Mouse IgM

Immunogen: Human articular cartilage aggrecan

Format:

100 μ L of ascites.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- ELISA
- Immunohistochemistry
- Immunoprecipitation
- Western blot

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.