

PRODUCT DATA SHEET

Product: Anti-Versican mAb, clone 5C12

Cat. No.: MC-1007 (100 µL)

Specificity:

Recognizes partially fucosylated 0-linked oligosaccharide epitopes prevailing on versican. Cross reactivity with aggrecans from different species.

Species Reactivity:

Human, bovine and chicken. Other species not tested.

Ig Isotype:

Mouse IgM

Immunogen:

Bovine aorta versican.

Format:

100 µL of ascites.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- Immunohistochemistry
- Western Blot
- ELISA

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.