

# PRODUCT DATA SHEET

Product: Anti-Versican mAb, clone 1H4

Cat. No.: MC-1008 (10 mL)

## Specificity:

Detects different isoforms of versican.

# Species Reactivity:

Human and bovine. Other species not tested.

# Ig Isotype:

Mouse IgG1

### Immunogen:

Bovine aorta versican.

#### Format:

10 mL of supernatant.

### Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

# Applications:

- Western Blot
- ELISA

The optimal dilution for a specific application should be determined by the researcher.

#### Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

### Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.