

PRODUCT DATA SHEET

Product: Anti-Versican mAb, clone 8D5

Cat. No.: MC-1009 (10 mL)

Specificity:

Directed against fucosyl capping structures of Nlinked oligosaccharide moieties uniquely expressed by vascular bovine versicans.

Species Reactivity:

Human and bovine. Other species not tested.

Ig Isotype:

Mouse IgG1

Immunogen:

Bovine aorta versican.

Format:

10 mL of supernatant.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- Immunohistochemistry
- Western Blot
- ELISA

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.