

PRODUCT DATA SHEET

Product: Anti-Laminin-5, -6 mAb, clone 667E12

Cat. No.: MC-1012 (100 µL)

Specificity:

Reacts with Laminin-5 and Laminin-6.

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Mouse IgG2a

Immunogen:

Recombinant human α 3 laminin chain fragment.

Format:

100 µL of ascites.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

 Immunohistochemistry (frozen sections): Use at a 1:50-1:100 dilution.

■ Western Blot: Use at a 1:500-1:1,000 dilution. The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.