



PRODUCT DATA SHEET

Product: Anti-Annexin V (VAC- α) mAb, clone VAA-33

Cat. No.: MC-250 (100 μ g)

Specificity:

Recognizes Annexin V (VAC- α)

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Mouse IgG_{2a/k}

Immunogen:

Recombinant Human Annexin V

Format:

100 μ g of purified monoclonal antibody at 1 mg/mL in PBS. >95% IgG_{2a/k} as determined by SDS-gel electrophoresis prior to addition of stabilizer.

Storage:

Store at 4°C. After opening, aliquot and freeze at -20°C.

Applications:

- ELISA
- Flow cytometry: Suitable for use as primary staining antibody. Use of MC-251, FITC conjugated antibody (sold separately), can be used for direct cell staining and is recommended for intracellular staining.
- Western blot

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.