



PRODUCT DATA SHEET

Product: Anti-Annexin A5 (Annexin V)-Biotin, Clone VAA-33

Cat. No: MC-252 (100 tests)

Specificity:

React with human Annexin V

Species Reactivity:

Human, others not tested.

Ig Isotype: Ig₂α/k

Immunogen:

Recombinant human Annexin V

Format:

Purified, biotin conjugated antibody in 1 ml of 50 mM Tris-HCl, 100 mM NaCl, 1% BSA, 0.02% sodium azide, pH 7.4 (>95% IgG2a as determined by SDS-gel electrophoresis prior to addition of stabilizer and preservative).

Storage and Stability:

Store at 2-8°C.

Applications and Suggested Dilutions:

- FACS analysis (Note: staining is intracellular, and cells must be fixed and permeabilized). This product may be used in a two-step procedure using the biotin-(strept-)avidin system.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.