

PRODUCT DATA SHEET

Product: Anti-CD40, Clone 3B2

Cat. No: MC-276 (100 μg)

Specificity:

This antibody reacts with CD40.

Species Reactivity:

Human, Others not tested.

Ig Isotype: Mouse IgG1/ κ

Format:

100 $\mu g/100$ μl of purified antibody in PBS.

Purified by protein A column.

Storage and Stability:

Store at -20 ℃.

Applications and Suggested Dilutions:

- ELISA: detection Ab, recommend 1 μg/ml.
- Immunohistochemistry: On cryostat sections, recommend 10 μg/ml.
- Western blot: recommend 10 μg/ml.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.