

PRODUCT DATA SHEET

Product: Anti-RANKL mAb, clone Ranky-1

Cat. No.: MC-357 (100 μg)

Synonyms:

ODF, OPGL, TRANCE, TNFSF11, CD254

Specificity:

Recognizes human RANKL.

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Rat IgG₁

Immunogen:

Recombinant human soluble RANKL extracellular domain (aa 151-316)

Format:

 $100~\mu g$ of antibody purified from concentrated hybridoma tissue culture supernatant. Antibody is at 1 mg/mL in PBS with 0.02% sodium azide.

Storage:

Store at 4°C for short term and at -20°C long term. Aliquot to avoid freeze/thaw cycles.

Applications and Suggested Dilutions:

- ELISA: Capture
- Immunohistochemistry: Paraffin sections
- Western blot: Use at a 1:1,000 dilution.

Note: Does **NOT** work for flow cytometry or immunoprecipitation.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.