



PRODUCT DATA SHEET

Product: Anti-SARM mAb, clone Sarmy-1

Cat. No.: MC-368 (100 µg)

Synonyms:

TIR-1, Sterile α and TIR Motif-containing Protein 1, SAM and ARM-containing Protein, Sterile α motif and Armadillo Repeat-containing protein.

Specificity:

Recognizes human SARM

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Mouse IgG_{2b}

Immunogen:

Recombinant human SARM (aa 93-292)

Format:

100 µg of antibody at 1 mg/mL purified from concentrated hybridoma tissue culture supernatant in PBS containing 10% glycerol with 0.02% sodium azide. $\geq 95\%$ purity by SDS-PAGE.

Storage:

Store at 4°C for short term and at -20°C long term. Aliquot to avoid freeze/thaw cycles.

Application and Suggested Dilutions:

■ Immunoprecipitation: Use at a 1:200 dilution.

Note: Does NOT work with Western blots.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.