



## PRODUCT DATA SHEET

**Product:** Anti-RIP3 mAb, clone Rippy-3

**Cat. No.:** MC-573 (100 µg)

**Synonyms:**

Receptor-interacting Serine/Threonine Protein Kinase 3; RIP-like Protein Kinase 3; Receptor-interacting Protein 3

**Specificity:**

Recognizes human RIP3.

**Species Reactivity:**

Human. Other species not tested.

**Ig Isotype:**

Mouse IgG<sub>1</sub>

**Immunogen:**

Recombinant human RIP3

**Format:**

100 µg of purified antibody at 1 mg/mL in PBS containing 0.02% sodium azide. Purity ≥95% as determined by SDS-PAGE.

**Storage and Stability:**

Store at 4°C short term and at -20°C long term. Aliquot to avoid freeze/thaw cycles.

**Application:**

■ Western blot

The researcher should determine the optimal dilution for a specific application.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.