



PRODUCT DATA SHEET

Product: Anti-Versican mAb, clone 2E12

Cat. No.: MC-809 (100 μ L)

Specificity:

Reacts with N-linked oligosaccharides specific of versicans. The antibody is not species-specific.

Ig Isotype:

Mouse IgG1

Immunogen:

Bovine vascular proteoglycans

Format:

100 μ L ascites fluid containing 0.01% sodium azide. Antibody concentration >1mg/mL.

Storage:

Store at -20°C long term, 4°C short term.

Applications and Suggested Dilutions:

- ELISA: 1:500 - 1:1,000
- Immunohistochemistry: (paraffin embedded): 1:50 - 1:3,000.
- Immunoprecipitation: 1:50 - 1:300
- Western Blotting: 1:100 - 1:300

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.