



PRODUCT DATA SHEET

Product: Anti-Laminin α 4 Subunit mAb, clone 652C4

Cat. No.: MC-818 (100 μ L)

Specificity:

Reacts with laminin α 4 subunit.

Species Reactivity:

Human. Others not tested.

Ig Isotype:

Mouse IgG₁.

Immunogen:

Recombinant fragment corresponding to the G1-3 domains of human laminin α 4 subunit.

Epitope:

Within the G1-3 domains.

Format:

100 μ L ascites fluid. Antibody concentration 1 - 3 mg/mL.

Storage:

Store at 4°C short term, -20°C long term. Aliquot to avoid cycles of freeze/thaw.

Applications and Suggested Dilutions:

- Immunohistochemistry:
(cryostat sections): 1:100 - 1:300.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.