

PRODUCT DATA SHEET

Product: Anti-Vimentin mAb, clone VI-RE/1

Cat. No.: MC-891 (100 µg)

Specificity: Reacts with Vimentin.

Ig Isotype:

Mouse IgG₁

Species Reactivity:

Human. Others not tested.

Format:

100 μg of monoclonal antibody at 1 mg/mL in PBS containing 15 mM sodium azide, pH 7.4.

Purity:

>95% by SDS-PAGE.

Storage:

Store at 4° C for short term. Aliquot and store at -20° C for long term. Avoid repeated freeze/thaw cycles.

Applications:

- ELISA
- Western blot
- Immunocytochemistry

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.