

PRODUCT DATA SHEET

Product: Anti-Keratin 4, clone 6B10

Cat. No: MC-959 (1 ml)

Specificity: Recognizes keratin 4 (59 kD)

Species Reactivity: Human, others not

tested

Isotype: Mouse IgG₁

Immunogen: Cytokeratin preparation purified

from human esophagus.

Format: 1 ml hybridoma sulture supernatant with approximately 15 ug antibody/ml and 1%

FCS and 10 mM sodium azide.

Storage and Stability: Store undiluted antibody at 2-8°C. For a long-term storage, aliquot and store at -20°C.

Applications and Suggested Dilutions:

- Immunohistochemistry –Frozen sections at 1:40 dilution, paraffin-embedded tissue sections at 1:10 dilution.
- Western blot
- Two dimensional gel electrophoresis followed by immunoblotting.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.