



PRODUCT DATA SHEET

Product: Anti-Creatine Kinase-MB mAb, clone 180

Cat. No: RM-096

Specificity:

Reacts with human creatine Kinase-MB

Species Reactivity:

Human. Others not tested.

Ig Isotype:

Mouse IgG₁

Purification:

Purified by Protein A chromatography

Format:

Purified antibody supplied in 15 mM potassium phosphate, 0.15 M NaCl, pH 7.2, with 0.09% sodium azide added as preservative

Storage and Stability:

Store antibody at 2 - 8°C short term. For long term, aliquot and store at -20°C. Avoid freeze/thaw cycles.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.