



PRODUCT DATA SHEET

Product: Anti-CD44, polyclonal

Cat. No.: PC-016 (100 µL)

Specificity:

Recognizes human CD44.

Species Reactivity:

Human. Other species not tested.

Host:

Rabbit

Immunogen:

Recombinant human CD44 from bacteria.

Format:

100 µL ascites fluid. Antibody concentration at 1-5 mg/mL.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- ELISA
- Immunocytochemistry
- Western Blots

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.