



## PRODUCT DATA SHEET

**Product:** Anti-MMP-13, polyclonal

**Cat. No.:** PC-130 (100  $\mu$ L)

**Specificity:**

Recognizes MMP-13. Does not cross-react with MMP-2, MMP-9, MMP-14, MMP-15, MMP-16 or MMP-17.

**Species Reactivity:**

Human, other species not tested.

**Ig Isotype:**

Rabbit IgG

**Immunogen:**

Human recombinant matrix metalloproteinase 13 (MMP-13).

**Format:**

100  $\mu$ L of purified antibody, contains 0.01% thimerosal.

**Storage and Stability:**

Store at 4°C.

**Applications:**

- ELISA
- Western blot

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.