

# PRODUCT DATA SHEET

**Product:** Anti-Beta Galactosidase pAb

Cat. No.: PC-215 (0.1 mg)

#### Procedure:

Rabbits were immunized with highly purified beta-galactosidase from *E. coli* and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

### Host:

Rabbit

### Format:

Liquid. 0.1 mL of affinity purified polyclonal antibody at 1.0 mg/mL in PBS, pH 7.2 with 0.09% sodium azide.

## Storage and Stability:

Store at 4°C. Stable for 24 months.

## Applications:

- ELISA
- Blotting
- IHC

The optimal dilution for a specific application should be determined by the researcher.

#### Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

## Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.