



## PRODUCT DATA SHEET

**Product:** Anti-C-MYC (EQKLISEEDL) pAb

**Cat. No.:** PC-221 (0.1 mg)

**Procedure:**

Rabbits were immunized with highly purified EQKLISEEDL conjugated with KLH (C-MYC) and the resulting antiserum was collected. Antibody was solid phase adsorbed then immunoaffinity purified off an antigen containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**Host:**

Rabbit

**Format:**

Liquid. 0.1 mL affinity purified antibody at 1 mg/mL in PBS, pH 7.2 with 0.09% sodium azide.

**Storage:**

Store at 4°C.

**Applications:**

- ELISA
- Immunohistochemistry
- Blotting

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.