



## PRODUCT DATA SHEET

**Product:** Anti-C-Myc (EQKLISEEDL), Polyclonal

**Cat. No:** PC-222 (1.0 mg)

**Procedure:**

Rabbits were immunized with the C-Myc peptide EQKLISEEDL conjugated with KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid phase.

**Host:** Rabbit

**Format:**

Liquid. Antibody in PBS, pH 7.2 with 0.1% sodium azide.

**Storage and Stability:**

Store at +4°C.

**Applications and Suggested Dilutions:**

- ELISA: coating 1:200 – 1:500
- Immunofluorescence: Use at 1:250 – 1:1,000
- Western blot: Use at 1:2,000-1:10,000.

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.