

PRODUCT DATA SHEET

Product: Anti-FLAG (DYKDDDDK) pAb

Cat. No.: PC-237 (0.1 mg)

Procedure:

Rabbits were immunized with the FLAG peptide conjugated to KLH. After multiple immunizations in Freund's adjuvant, serum was collected and the antibodies were affinity-purified using the peptide coupled to a different carrier.

Host:

Rabbit

Format:

Liquid. 100 μ L of affinity-purified polyclonal antibody at 1.0 mg/mL in PBS, pH 7.2, with 0.1% sodium azide.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA: For detection use at 0.04-0.2 μg/mL. For coating start at 5.0 μg/mL.
- Western blot: Use at 0.2-1.0 μg/mL.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.