



PRODUCT DATA SHEET

Product: Anti-Hemagglutinin (HA), FITC-labeled, Polyclonal

Cat. No: PC-261 (1.0 mg)

Procedure:

Rabbits were immunized with HA-tag conjugated to KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide conjugated to a solid-phase support and the antibodies were conjugated with fluorescein isothiocyanate (FITC).

Host: Rabbit

Specification:

Molar F/P ratio = 3.7.

Format:

Liquid. FITC-labeled antibody in PBS, pH 7.2 with 0.2% BSA and 0.1% sodium azide.

Storage and Stability:

Store at +4°C.

Applications and Suggested Dilutions:

- ELISA: 1:5,000 – 1:15,000
- Immunofluorescence: 1:250 – 1:1,000
- Western blot: 1:2,000 – 1:10,000

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.