



PRODUCT DATA SHEET

Product: Anti-Hexahistidine pAb

Cat. No.: PC-275 (0.1 mg)

Procedure:

Rabbits were immunized with highly purified HHHHHH-KLH (6xHIS) and the resulting antiserum was collected. Antibody was solid phase absorbed and then immunoaffinity purified off an antigen containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

Host:

Rabbit

Format:

Liquid. 0.1 mL of affinity-purified antibody at 1.0 mg/mL in PBS, pH 7.2 and 0.09% sodium azide.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA
- Blotting
- IHC

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product.

KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.