



PRODUCT DATA SHEET

Product: Anti-V5 (GKPIPPELLGLDST) pAb

Cat. No.: PC-300 (0.1 mg)

Procedure:

Rabbits were immunized with highly purified V5 peptide conjugated with KLH and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG

Host:

Rabbit

Format:

Liquid. 100 μ L of affinity-purified polyclonal antibody at 1 mg/mL in PBS, pH 7.2 with 0.09% sodium azide.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA
- Immunohistochemistry
- Blotting

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.