



PRODUCT DATA SHEET

Product: Anti-C-MYC (EQKLISEEDL)-FITC pAb

Cat. No.: PC-308 (0.1 mg)

Procedure:

Chickens were immunized with C-MYC (EQKLISEEDL) peptide conjugated with KLH. After multiple immunizations in Freund's adjuvant, eggs were collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid phase. The antibodies were conjugated with Fluorescein Isothiocyanate (FITC).

Specification:

Prior to conjugation, the antibodies were analyzed by ELISA and Western blot analyses. The molar F/P ratio is 5.6.

Host:

Chicken

Format:

100 μ L of 1.0 mg/mL liquid FITC-labeled polyclonal antibody in PBS, pH 7.2 with protein stabilizer and 0.1% sodium azide.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA
- Western blot: Use at a 1:1,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.