

PRODUCT DATA SHEET

Product: Anti-C-MYC-HRP pAb

Cat. No.: PC-309 (0.1 mg)

Procedure:

Chickens were immunized with the C-MYC peptide conjugated with KLH. After multiple immunizations in Freund's adjuvant, eggs were collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid phase. The antibodies were conjugated with horseradish peroxidase (HRP).

Host:

Chicken

Format:

Liquid. 0.1 mL of HRP-conjugated antibody at 1.0 mg/mL in 50 mM HEPES, saline with stabilizer and 0.002% thimerosal.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA: Use at 0.2 1.0 μg/mL.
- Western blot: The antibody was an effective probe at 0.5 μg/mL.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.