



PRODUCT DATA SHEET

Product: Anti-C-Myc (EQKLISEEDL)-HRP pAb

Cat. No.: PC-311 (0.1 mg)

Procedure:

Rabbits were immunized with the C-Myc peptide conjugated with KLH. After multiple immunizations, serum was collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid phase and conjugated with horseradish peroxidase (HRP).

Host:

Rabbit

Format:

Liquid. HRP-conjugated antibody at 1.0 mg/mL in 50 mM HEPES - saline with stabilizer and 0.002% thimerosal.

Storage and Stability:

Store at 4°C. Stable for 12 months after shipment date.

Applications and Suggested Dilutions:

- ELISA: Use at 1:5,000 – 1:20,000 for coating.
- Western blot: Use at 0.2 - 1.0 µg/mL.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.