



PRODUCT DATA SHEET

Product: Anti-S-tag pAb

Cat. No.: PC-335 (0.1 mg)

Procedure:

Rabbits were immunized with highly purified S-tag peptide conjugated with KLHA and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

Host:

Rabbit

Format:

Liquid. 0.1 mg of affinity purified polyclonal antibody at 1.0 mg/mL in PBS, pH 7.2 with 0.09% sodium azide.

Storage:

Store at 4°C. Stable for 24 months.

Applications and Suggested Dilutions:

- ELISA
- Blotting
- Immunohistochemistry

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.