



PRODUCT DATA SHEET

Product: Anti-S-tag, Polyclonal

Cat. No: PC-336 (1.0 mg)

Procedure:

Rabbits were immunized with S-tag conjugated with KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using peptide conjugated with a different carrier.

Host: Rabbit

Format:

Liquid. Purified antibody in PBS, pH 7.2 and 0.1% sodium azide.

Storage and Stability:

Store at +4°C.

Applications and Suggested Dilutions:

- ELISA: Use at 1:200 – 1:500 (coating);
1:1,000 – 1:20,000 (detection)
- Western blot: Use at 1:1,000 – 1:5,000

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.