



PRODUCT DATA SHEET

Product: Anti-S-tag-HRP pAb

Cat. No.: PC-340 (0.1 mg)

Procedure:

Rabbits were immunized with highly purified S-tag peptide conjugated with KLH and the resulting antiserum was collected. Antibody was solid phase adsorbed, immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

Host:

Rabbit

Format:

Liquid. 0.1 mL of HRP-labeled antibody at 1 mg/mL in HEPES - Saline with stabilizer and 0.002% thimerosal.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA
- Blotting
- Immunohistochemistry

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.