

# **PRODUCT DATA SHEET**

## Product: Anti-Vitronectin, polyclonal

Cat. No.: PC-398 (100 µL)

*Specificity:* Recognizes human vitronectin.

#### Species Reactivity:

Human. Other species not tested.

*Host:* Rabbit

*Immunogen:* Human vitronectin.

*Format:* 100 μL of rabbit antiserum.

#### Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

### Applications:

- ELISA
- Immunocytochemistry
- Western Blots

The optimal dilution for a specific application should be determined by the researcher.

#### Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

#### Warranty:

No warranties, expressed or implied, are made regarding the use of this product. **KAMIYA BIOMEDICAL COMPANY** is not liable for any damage, personal injury, or economic loss caused by this product.