

# PRODUCT DATA SHEET

**Product:** Anti-PYY pAb

Cat. No.: PC-509 (50 µL Lyophilized)

## Immunogen:

Synthetic pig peptide YY (PYY).

# Species Reactivity:

Mouse, rat and cat. Other species not tested.

#### Host:

Guinea Pig

## Positive Controls:

Frozen sections of rat intestine. Absorption with 10-100 µg PYY per mL diluted pAb abolishes the staining.

#### Format:

50 µL undiluted, lyophilized polyclonal antibody.

## Reconstitution:

Dissolve the lyophilized antiserum in 50-100  $\mu$ L distilled water, and dilute further in 0.1 M PBS with 1% BSA and 0.1% NaN<sub>3</sub>.

#### Storage:

Store lyophilized pAb at 4°C until expiration date. Aliquot and store reconstituted pAb at -20°C. Avoid repeated freeze/thaw cycles.

# Applications and Suggested Dilution:

■ Immunohistochemistry: Paraffin and frozen sections. Use at a 1:1,500 – 1:2,500 dilution using FITC, with over-night incubation at 4°C.

The optimal dilution for a specific application should be determined by the researcher.

#### Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

#### Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.