



PRODUCT DATA SHEET

Product: Anti-Laminin-1 pAb

Cat. No.: PC-562 (100 μ L)

Specificity:

Reacts with the laminin α 1 chain and weakly reacts with β 1 and γ 1 chains.

Species Reactivity:

Human, bovine, mouse and rat. Other species not tested.

Ig Isotype:

Rabbit IgG

Immunogen:

Purified laminin-1 from EHS tumor.

Format:

100 μ L of antiserum.

Storage:

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- Western Blot
- ELISA

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.