



## PRODUCT DATA SHEET

**Product:** Anti-Glu-Glu (EYMPME) pAb

**Cat. No.:** PC-246 (0.1 mg)

**Procedure:**

Rabbits were immunized with the Glu-Glu (EYMPME) conjugated with KLH. After multiple immunizations in Freund's adjuvant, the serum was collected and the antibodies were immunoaffinity-purified using the peptide conjugated with a different carrier.

**Host:**

Rabbit

**Format:**

Liquid. Affinity-purified antibody at 1.0 mg/mL in PBS, pH 7.2, with 0.1% sodium azide.

**Storage:**

Store at 4°C.

**Applications and Suggested Dilutions:**

- ELISA (Indirect): Start at 1.0 µg/mL.
- Blotting: Use at 0.2 - 1 µg/mL.

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.