

PRODUCT DATA SHEET

Product: Anti-Hemagglutinin (HA) pAb

Cat. No.: PC-260 (1.0 mg)

Procedure:

Rabbits were immunized with YPYDVPDYA (influenza hemagglutinin-HA-epitope) conjugated to KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide conjugated to a solid-phase support.

Host:

Rabbit

Format:

Liquid. 1.0 mL of affinity-purified polyclonal antibody at 1.0 mg/mL in PBS, pH 7.2, and 0.1% sodium azide.


Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA: Use at a 1:200 – 1:500 dilution (capture); 1:5,000 - 1:20,000 dilution (detection).
- Immunofluorescence: Use at a 1:250 – 1:2,000 dilution.
- Western blot: Use at a 1:5,000 – 1:10,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.



Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.